





# LEFTY1/LEFTY2 Antibody

<b>Product Code</b>	CSB-PA006317
<b>Storage</b>	Upon receipt, store at -20°C or -80°C. Avoid repeated freeze.
<b>Uniprot No.</b>	O75610/O00292
<b>Immunogen</b>	Synthesized peptide derived from the C-terminal region of Human Lefty.
<b>Raised In</b>	Rabbit
<b>Species Reactivity</b>	Human, Mouse
<b>Tested Applications</b>	WB, ELISA; WB:1:500-1:2000, ELISA:1:10000
<b>Form</b>	Liquid
<b>Conjugate</b>	Non-conjugated
<b>Storage Buffer</b>	Liquid in PBS containing 50% glycerol, 0.5% BSA and 0.02% sodium azide.
<b>Purification Method</b>	The antibody was affinity-purified from rabbit antiserum by affinity-chromatography using epitope-specific immunogen.
<b>Isotype</b>	IgG
<b>Alias</b>	LEFTY1; LEFTB; LEFTYB; Left-right determination factor 1; Left-right determination factor B; Protein lefty-1; Protein lefty-B; LEFTY2; EBAF; LEFTA; LEFTYA; TGFB4; PSEC0024; Left-right determination factor 2; Endometrial bleeding-associated factor; Left-ri
<b>Immunogen Species</b>	Homo sapiens (Human)
<b>Target Names</b>	LEFTY1/LEFTY2

## Image



Western Blot analysis of HEPG2-UV cells using Lefty Polyclonal Antibody



Western Blot analysis of HepG2-UV cells using Lefty Polyclonal Antibody

## Usage

For Research Use Only. Not for use in diagnostic or therapeutic procedures.