





LEFTY1/LEFTY2 Antibody

Product Code	CSB-PA006317
Storage	Upon receipt, store at -20°C or -80°C. Avoid repeated freeze.
Uniprot No.	O75610/O00292
Immunogen	Synthesized peptide derived from the C-terminal region of Human Lefty.
Raised In	Rabbit
Species Reactivity	Human, Mouse
Tested Applications	WB, ELISA; WB:1:500-1:2000, ELISA:1:10000
Form	Liquid
Conjugate	Non-conjugated
Storage Buffer	Liquid in PBS containing 50% glycerol, 0.5% BSA and 0.02% sodium azide.
Purification Method	The antibody was affinity-purified from rabbit antiserum by affinity-chromatography using epitope-specific immunogen.
Isotype	IgG
Alias	LEFTY1; LEFTB; LEFTYB; Left-right determination factor 1; Left-right determination factor B; Protein lefty-1; Protein lefty-B; LEFTY2; EBAF; LEFTA; LEFTYA; TGFB4; PSEC0024; Left-right determination factor 2; Endometrial bleeding-associated factor; Left-ri
Immunogen Species	Homo sapiens (Human)
Target Names	LEFTY1/LEFTY2

Image



Western Blot analysis of HEPG2-UV cells using Lefty Polyclonal Antibody



Western Blot analysis of HepG2-UV cells using Lefty Polyclonal Antibody