





NRTN Antibody

Product Code	CSB-PA006329
Storage	Upon receipt, store at -20°C or -80°C. Avoid repeated freeze.
Uniprot No.	Q99748
Immunogen	Synthesized peptide derived from the Internal region of Human NTN.
Raised In	Rabbit
Species Reactivity	Human
Tested Applications	WB, ELISA; WB:1:500-1:2000, ELISA:1:10000
Form	Liquid
Conjugate	Non-conjugated
Storage Buffer	Liquid in PBS containing 50% glycerol, 0.5% BSA and 0.02% sodium azide.
Purification Method	The antibody was affinity-purified from rabbit antiserum by affinity-chromatography using epitope-specific immunogen.
Isotype	IgG
Alias	NRTN; Neurturin
Immunogen Species	Homo sapiens (Human)
Target Names	NRTN

Image



Western Blot analysis of K562 cells using NTN Polyclonal Antibody



Western Blot analysis of K562 L929 cells using NTN Polyclonal Antibody



Usage

For Research Use Only. Not for use in diagnostic or therapeutic procedures.