





EPGN Antibody

Product Code	CSB-PA006336
Storage	Upon receipt, store at -20°C or -80°C. Avoid repeated freeze.
Uniprot No.	Q6UW88
Immunogen	Synthesized peptide derived from the Internal region of Human Epigen.
Raised In	Rabbit
Species Reactivity	Human
Tested Applications	WB, ELISA; WB:1:500-1:2000, ELISA:1:10000
Form	Liquid
Conjugate	Non-conjugated
Storage Buffer	Liquid in PBS containing 50% glycerol, 0.5% BSA and 0.02% sodium azide.
Purification Method	The antibody was affinity-purified from rabbit antiserum by affinity-chromatography using epitope-specific immunogen.
Isotype	IgG
Alias	EPGN; Epigen; Epithelial mitogen; EPG
Immunogen Species	Homo sapiens (Human)
Target Names	EPGN

Image



Western Blot analysis of 293T cells using Epigen Polyclonal Antibody



Western Blot analysis of 293T 823 cells using Epigen Polyclonal Antibody



Usage

For Research Use Only. Not for use in diagnostic or therapeutic procedures.