





FUNDC1 Antibody

Product Code	CSB-PA006338
Storage	Upon receipt, store at -20°C or -80°C. Avoid repeated freeze.
Uniprot No.	Q8IVP5
Immunogen	Synthesized peptide derived from the N-terminal region of Human FUNDC1.
Raised In	Rabbit
Species Reactivity	Human, Mouse, Rat
Tested Applications	WB, ELISA; WB:1:500-1:2000, ELISA:1:10000
Form	Liquid
Conjugate	Non-conjugated
Storage Buffer	Liquid in PBS containing 50% glycerol, 0.5% BSA and 0.02% sodium azide.
Purification Method	The antibody was affinity-purified from rabbit antiserum by affinity-chromatography using epitope-specific immunogen.
Isotype	IgG
Alias	FUNDC1; FUN14 domain-containing protein 1
Immunogen Species	Homo sapiens (Human)
Target Names	FUNDC1

Image



Western Blot analysis of AD293 cells using FUNDC1 Polyclonal Antibody



Western Blot analysis of AD293 cells using FUNDC1 Polyclonal Antibody



Usage

For Research Use Only. Not for use in diagnostic or therapeutic procedures.