




Phospho-ATP1A1 (S23) Antibody

Product Code	CSB-PA006343
Storage	Upon receipt, store at -20°C or -80°C. Avoid repeated freeze.
Uniprot No.	P05023
Immunogen	Synthesized peptide derived from Human Na ⁺ /K ⁺ -ATPase α1 around the phosphorylation site of S23.
Raised In	Rabbit
Species Reactivity	Rat
Tested Applications	WB, IF, ELISA; WB:1:500-1:2000, IF:1:200-1:1000, ELISA:1:20000
Form	Liquid
Conjugate	Non-conjugated
Storage Buffer	Liquid in PBS containing 50% glycerol, 0.5% BSA and 0.02% sodium azide.
Purification Method	The antibody was affinity-purified from rabbit antiserum by affinity-chromatography using epitope-specific immunogen.
Isotype	IgG
Alias	ATP1A1; Sodium/potassium-transporting ATPase subunit alpha-1; Na(+)/K(+) ATPase alpha-1 subunit; Sodium pump subunit alpha-1
Immunogen Species	Rattus norvegicus (Rat)
Target Names	ATP1A1

Image



Western Blot analysis of rat cells using Phospho-Na⁺/K⁺-ATPase α1 (S23) Polyclonal Antibody

Usage

For Research Use Only. Not for use in diagnostic or therapeutic procedures.