



RASGRF1 Antibody

| | |
|----------------------------|---|
| Product Code | CSB-PA006353 |
| Storage | Upon receipt, store at -20°C or -80°C. Avoid repeated freeze. |
| Uniprot No. | Q13972 |
| Immunogen | Synthesized peptide derived from Mouse Ras-GRF1 around the non-phosphorylation site of S916. |
| Raised In | Rabbit |
| Species Reactivity | Mouse,Rat |
| Tested Applications | IHC, ELISA; IHC:1:100-1:300, ELISA:1:5000 |
| Form | Liquid |
| Conjugate | Non-conjugated |
| Storage Buffer | Liquid in PBS containing 50% glycerol, 0.5% BSA and 0.02% sodium azide. |
| Purification Method | The antibody was affinity-purified from rabbit antiserum by affinity-chromatography using epitope-specific immunogen. |
| Isotype | IgG |
| Alias | RASGRF1; CDC25; GNRP; GRF1; Ras-specific guanine nucleotide-releasing factor 1; Ras-GRF1; Guanine nucleotide-releasing protein; GNRP; Ras-specific nucleotide exchange factor CDC25 |
| Immunogen Species | Mus musculus (Mouse) |
| Target Names | RASGRF1 |
| Usage | For Research Use Only. Not for use in diagnostic or therapeutic procedures. |