




TNNI3 Antibody

Product Code	CSB-PA006356
Storage	Upon receipt, store at -20°C or -80°C. Avoid repeated freeze.
Uniprot No.	P19429
Immunogen	Synthesized peptide derived from Mouse Troponin I-C around the non-phosphorylation site of S22/S23.
Raised In	Rabbit
Species Reactivity	Mouse,Rat
Tested Applications	WB, IHC, IF, ELISA; WB:1:500-1:2000, IHC:1:100-1:300, IF:1:200-1:1000, ELISA:1:5000
Form	Liquid
Conjugate	Non-conjugated
Storage Buffer	Liquid in PBS containing 50% glycerol, 0.5% BSA and 0.02% sodium azide.
Purification Method	The antibody was affinity-purified from rabbit antiserum by affinity-chromatography using epitope-specific immunogen.
Isotype	IgG
Alias	TNNI3; TNNC1; Troponin I; cardiac muscle; Cardiac troponin I
Immunogen Species	Mus musculus (Mouse)
Target Names	TNNI3

Image



Western Blot analysis of mouse cells using Troponin I-C Polyclonal Antibody

Usage

For Research Use Only. Not for use in diagnostic or therapeutic procedures.