




Phospho-HDAC6 (S22) Antibody

Product Code	CSB-PA006360
Storage	Upon receipt, store at -20°C or -80°C. Avoid repeated freeze.
Uniprot No.	Q9UBN7
Immunogen	Synthesized peptide derived from Human HDAC6 around the phosphorylation site of S22.
Raised In	Rabbit
Species Reactivity	Human, Mouse
Tested Applications	WB, IHC, IF, ELISA; WB:1:500-1:2000, IHC:1:100-1:300, IF:1:200-1:1000, ELISA:1:10000
Form	Liquid
Conjugate	Non-conjugated
Storage Buffer	Liquid in PBS containing 50% glycerol, 0.5% BSA and 0.02% sodium azide.
Purification Method	The antibody was affinity-purified from rabbit antiserum by affinity-chromatography using epitope-specific immunogen.
Isotype	IgG
Alias	HDAC6; KIAA0901; JM21; Histone deacetylase 6; HD6
Immunogen Species	Homo sapiens (Human)
Target Names	HDAC6

Image



Western Blot analysis of 3T3 cells using Phospho-HDAC6 (S22) Polyclonal Antibody