





CDH5 Antibody

Product Code	CSB-PA006367
Storage	Upon receipt, store at -20°C or -80°C. Avoid repeated freeze.
Uniprot No.	P33151
Immunogen	Synthesized peptide derived from Human VE-Cadherin around the non-phosphorylation site of Y731.
Raised In	Rabbit
Species Reactivity	Human, Mouse
Tested Applications	WB, IHC, ELISA; WB:1:500-1:2000, IHC:1:100-1:300, ELISA:1:40000
Form	Liquid
Conjugate	Non-conjugated
Storage Buffer	Liquid in PBS containing 50% glycerol, 0.5% BSA and 0.02% sodium azide.
Purification Method	The antibody was affinity-purified from rabbit antiserum by affinity-chromatography using epitope-specific immunogen.
Isotype	IgG
Alias	CDH5; Cadherin-5; 7B4 antigen; Vascular endothelial cadherin; VE-cadherin; CD antigen CD144
Immunogen Species	Homo sapiens (Human)
Target Names	CDH5

Image



Western Blot analysis of Jurkat cells using VE-Cadherin Polyclonal Antibody



Western Blot analysis of A549 using VE-Cadherin Polyclonal Antibody.

Usage

For Research Use Only. Not for use in diagnostic or therapeutic procedures.