




Phospho-DNM1L (S637) Antibody

Product Code	CSB-PA006377
Storage	Upon receipt, store at -20°C or -80°C. Avoid repeated freeze.
Uniprot No.	O00429
Immunogen	Synthesized peptide derived from Human DRP1 around the phosphorylation site of S637.
Raised In	Rabbit
Species Reactivity	Human,Mouse,Rat
Tested Applications	WB, IHC, ELISA; WB:1:500-1:2000, IHC:1:100-1:300, ELISA:1:10000
Form	Liquid
Conjugate	Non-conjugated
Storage Buffer	Liquid in PBS containing 50% glycerol, 0.5% BSA and 0.02% sodium azide.
Purification Method	The antibody was affinity-purified from rabbit antiserum by affinity-chromatography using epitope-specific immunogen.
Isotype	IgG
Alias	DNM1L; DLP1; DRP1; Dynamin-1-like protein; Dnm1p/Vps1p-like protein; DVLP; Dynamin family member proline-rich carboxyl-terminal domain less; Dymple; Dynamin-like protein; Dynamin-like protein 4; Dynamin-like protein IV; HdynIV; Dynamin-rela
Immunogen Species	Homo sapiens (Human)
Target Names	DNM1L

Image



Western Blot analysis of JK cells using Phospho-DRP1 (S637) Polyclonal Antibody

Usage

For Research Use Only. Not for use in diagnostic or therapeutic procedures.