



# DNAL4 Antibody

|                            |   |
|----------------------------|---|
| <b>Product Code</b>        | CSB-PA006388  |
| <b>Storage</b>             | Upon receipt, store at -20°C or -80°C. Avoid repeated freeze.   |
| <b>Uniprot No.</b>         | O96015  |
| <b>Immunogen</b>           | Synthesized peptide derived from the N-terminal region of Human Dynein LC 4.  |
| <b>Raised In</b>           | Rabbit  |
| <b>Species Reactivity</b>  | Human, Mouse  |
| <b>Tested Applications</b> | IHC, ELISA; IHC:1:100-1:300, ELISA:1:40000  |
| <b>Form</b>                | Liquid  |
| <b>Conjugate</b>           | Non-conjugated  |
| <b>Storage Buffer</b>      | Liquid in PBS containing 50% glycerol, 0.5% BSA and 0.02% sodium azide.   |
| <b>Purification Method</b> | The antibody was affinity-purified from rabbit antiserum by affinity-chromatography using epitope-specific immunogen. |
| <b>Isotype</b>             | IgG   |
| <b>Alias</b>               | DNAL4; Dynein light chain 4; axonemal   |
| <b>Immunogen Species</b>   | Homo sapiens (Human)  |
| <b>Target Names</b>        | DNAL4   |