




Phospho-TOB1 (S164) Antibody

| | |
|----------------------------|---|
| Product Code | CSB-PA006393 |
| Storage | Upon receipt, store at -20°C or -80°C. Avoid repeated freeze. |
| Uniprot No. | P50616 |
| Immunogen | Synthesized peptide derived from Human TOB1 around the phosphorylation site of S164. |
| Raised In | Rabbit |
| Species Reactivity | Human,Mouse,Rat |
| Tested Applications | WB, IHC, ELISA; WB:1:500-1:2000, IHC:1:100-1:300, ELISA:1:5000 |
| Form | Liquid |
| Conjugate | Non-conjugated |
| Storage Buffer | Liquid in PBS containing 50% glycerol, 0.5% BSA and 0.02% sodium azide. |
| Purification Method | The antibody was affinity-purified from rabbit antiserum by affinity-chromatography using epitope-specific immunogen. |
| Isotype | IgG |
| Alias | TOB1; TOB; TROB1; Protein Tob1; Transducer of erbB-2 1 |
| Immunogen Species | Homo sapiens (Human) |
| Target Names | TOB1 |

Image



Western Blot analysis of JK cells using Phospho-TOB1 (S164) Polyclonal Antibody