




TNK2 Antibody

Product Code	CSB-PA006415
Storage	Upon receipt, store at -20°C or -80°C. Avoid repeated freeze.
Uniprot No.	Q07912
Immunogen	Synthesized peptide derived from Human ACK around the non-phosphorylation site of Y284.
Raised In	Rabbit
Species Reactivity	Human, Mouse
Tested Applications	WB, ELISA; WB:1:500-1:2000, ELISA:1:10000
Form	Liquid
Conjugate	Non-conjugated
Storage Buffer	Liquid in PBS containing 50% glycerol, 0.5% BSA and 0.02% sodium azide.
Purification Method	The antibody was affinity-purified from rabbit antiserum by affinity-chromatography using epitope-specific immunogen.
Isotype	IgG
Alias	TNK2; ACK1; Activated CDC42 kinase 1; ACK-1; Tyrosine kinase non-receptor protein 2
Immunogen Species	Homo sapiens (Human)
Target Names	TNK2

Image



Western Blot analysis of A549 cells using ACK Polyclonal Antibody

Usage

For Research Use Only. Not for use in diagnostic or therapeutic procedures.