




CDK5 Antibody

Product Code	CSB-PA006427
Storage	Upon receipt, store at -20°C or -80°C. Avoid repeated freeze.
Uniprot No.	Q00535
Immunogen	Synthesized peptide derived from Human Cdk5 around the non-phosphorylation site of Y15.
Raised In	Rabbit
Species Reactivity	Human,Mouse,Rat
Tested Applications	WB, IHC, ELISA; WB:1:500-1:2000, IHC:1:100-1:300, ELISA:1:5000
Form	Liquid
Conjugate	Non-conjugated
Storage Buffer	Liquid in PBS containing 50% glycerol, 0.5% BSA and 0.02% sodium azide.
Purification Method	The antibody was affinity-purified from rabbit antiserum by affinity-chromatography using epitope-specific immunogen.
Isotype	IgG
Alias	CDK5; CDKN5; Cyclin-dependent kinase 5; Cell division protein kinase 5; Serine/threonine-protein kinase PSSALRE; Tau protein kinase II catalytic subunit; TPKII catalytic subunit
Immunogen Species	Homo sapiens (Human)
Target Names	CDK5

Image



Western Blot analysis of HepG2 cells using Cdk5 Polyclonal Antibody

Usage For Research Use Only. Not for use in diagnostic or therapeutic procedures.