



ADAR Antibody

Product Code	CSB-PA006478
Storage	Upon receipt, store at -20°C or -80°C. Avoid repeated freeze.
Uniprot No.	P55265
Immunogen	Synthesized peptide derived from the C-terminal region of Human ADAR1.
Raised In	Rabbit
Species Reactivity	Human, Mouse, Rat
Tested Applications	IHC, ELISA; IHC:1:100-1:300, ELISA:1:20000
Form	Liquid
Conjugate	Non-conjugated
Storage Buffer	Liquid in PBS containing 50% glycerol, 0.5% BSA and 0.02% sodium azide.
Purification Method	The antibody was affinity-purified from rabbit antiserum by affinity-chromatography using epitope-specific immunogen.
Isotype	IgG
Alias	ADAR; ADAR1; DSRAD; G1P1; IFI4; Double-stranded RNA-specific adenosine deaminase; DRADA; 136 kDa double-stranded RNA-binding protein; p136; Interferon-inducible protein 4; IFI-4; K88DSRBP
Immunogen Species	Homo sapiens (Human)
Target Names	ADAR
Usage	For Research Use Only. Not for use in diagnostic or therapeutic procedures.