



SNAPC5 Antibody

Product Code	CSB-PA006481
Storage	Upon receipt, store at -20°C or -80°C. Avoid repeated freeze.
Uniprot No.	O75971
Immunogen	Synthesized peptide derived from the N-terminal region of Human SNAPC 19.
Raised In	Rabbit
Species Reactivity	Human, Mouse
Tested Applications	IHC, ELISA; IHC:1:100-1:300, ELISA:1:10000
Form	Liquid
Conjugate	Non-conjugated
Storage Buffer	Liquid in PBS containing 50% glycerol, 0.5% BSA and 0.02% sodium azide.
Purification Method	The antibody was affinity-purified from rabbit antiserum by affinity-chromatography using epitope-specific immunogen.
Isotype	IgG
Alias	SNAPC5; SNAP19; snRNA-activating protein complex subunit 5; SNAPc subunit 5; Small nuclear RNA-activating complex polypeptide 5; snRNA-activating protein complex 19 kDa subunit; SNAPc 19 kDa subunit
Immunogen Species	Homo sapiens (Human)
Target Names	SNAPC5
Usage	For Research Use Only. Not for use in diagnostic or therapeutic procedures.