




NDC80 Antibody

Product Code	CSB-PA006498
Storage	Upon receipt, store at -20°C or -80°C. Avoid repeated freeze.
Uniprot No.	O14777
Immunogen	Synthesized peptide derived from the Internal region of Human Hec1.
Raised In	Rabbit
Species Reactivity	Human, Mouse
Tested Applications	WB, IF, ELISA; WB:1:500-1:2000, IF:1:200-1:1000, ELISA:1:20000
Form	Liquid
Conjugate	Non-conjugated
Storage Buffer	Liquid in PBS containing 50% glycerol, 0.5% BSA and 0.02% sodium azide.
Purification Method	The antibody was affinity-purified from rabbit antiserum by affinity-chromatography using epitope-specific immunogen.
Isotype	IgG
Alias	NDC80; HEC; HEC1; KNTC2; Kinetochores protein NDC80 homolog; Highly expressed in cancer protein; Kinetochores protein Hec1; HsHec1; Kinetochores-associated protein 2; Retinoblastoma-associated protein HEC
Immunogen Species	Homo sapiens (Human)
Target Names	NDC80

Image



Western Blot analysis of Jurkat cells using Hec1 Polyclonal Antibody

Usage

For Research Use Only. Not for use in diagnostic or therapeutic procedures.