




SPON2 Antibody

Product Code	CSB-PA006509
Storage	Upon receipt, store at -20°C or -80°C. Avoid repeated freeze.
Uniprot No.	Q9BUD6
Immunogen	Synthesized peptide derived from the Internal region of Human Mindin.
Raised In	Rabbit
Species Reactivity	Human, Mouse, Rat, Monkey
Tested Applications	WB, ELISA; WB:1:500-1:2000, ELISA:1:40000
Form	Liquid
Conjugate	Non-conjugated
Storage Buffer	Liquid in PBS containing 50% glycerol, 0.5% BSA and 0.02% sodium azide.
Purification Method	The antibody was affinity-purified from rabbit antiserum by affinity-chromatography using epitope-specific immunogen.
Isotype	IgG
Alias	SPON2; DIL1; Spondin-2; Differentially expressed in cancerous and non-cancerous lung cells 1; DIL-1; Mindin
Immunogen Species	Homo sapiens (Human)
Target Names	SPON2

Image



Western Blot analysis of COS-7 cells using Mindin Polyclonal Antibody

Usage

For Research Use Only. Not for use in diagnostic or therapeutic procedures.