




N4BP2L2 Antibody

Product Code	CSB-PA006517
Storage	Upon receipt, store at -20°C or -80°C. Avoid repeated freeze.
Uniprot No.	Q92802
Immunogen	Synthesized peptide derived from the Internal region of Human PFAAP5.
Raised In	Rabbit
Species Reactivity	Human,Mouse
Tested Applications	WB, ELISA; WB:1:500-1:2000, ELISA:1:10000
Form	Liquid
Conjugate	Non-conjugated
Storage Buffer	Liquid in PBS containing 50% glycerol, 0.5% BSA and 0.02% sodium azide.
Purification Method	The antibody was affinity-purified from rabbit antiserum by affinity-chromatography using epitope-specific immunogen.
Isotype	IgG
Alias	N4BP2L2; CG005; PFAAP5; NEDD4-binding protein 2-like 2; Phosphonoformate immuno-associated protein 5
Immunogen Species	Homo sapiens (Human)
Target Names	N4BP2L2

Image



Western Blot analysis of LOVO cells using PFAAP5 Polyclonal Antibody