




Phospho-VAV3 (Y173) Antibody

Product Code	CSB-PA006523
Storage	Upon receipt, store at -20°C or -80°C. Avoid repeated freeze.
Uniprot No.	Q9UKW4
Immunogen	Synthesized peptide derived from Human Vav3 around the phosphorylation site of Y173.
Raised In	Rabbit
Species Reactivity	Human, Mouse
Tested Applications	WB, IHC, ELISA; WB:1:500-1:2000, IHC:1:100-1:300, ELISA:1:40000
Form	Liquid
Conjugate	Non-conjugated
Storage Buffer	Liquid in PBS containing 50% glycerol, 0.5% BSA and 0.02% sodium azide.
Purification Method	The antibody was affinity-purified from rabbit antiserum by affinity-chromatography using epitope-specific immunogen.
Isotype	IgG
Alias	VAV3; Guanine nucleotide exchange factor VAV3; VAV-3
Immunogen Species	Homo sapiens (Human)
Target Names	VAV3

Image



Western Blot analysis of JK cells using Phospho-Vav3 (Y173) Polyclonal Antibody

Usage

For Research Use Only. Not for use in diagnostic or therapeutic procedures.