




VAV3 Antibody

Product Code	CSB-PA006526
Storage	Upon receipt, store at -20°C or -80°C. Avoid repeated freeze.
Uniprot No.	Q9UKW4
Immunogen	Synthesized peptide derived from Human Vav3 around the non-phosphorylation site of Y173.
Raised In	Rabbit
Species Reactivity	Human, Mouse
Tested Applications	WB, ELISA; WB:1:500-1:2000, ELISA:1:5000
Form	Liquid
Conjugate	Non-conjugated
Storage Buffer	Liquid in PBS containing 50% glycerol, 0.5% BSA and 0.02% sodium azide.
Purification Method	The antibody was affinity-purified from rabbit antiserum by affinity-chromatography using epitope-specific immunogen.
Isotype	IgG
Alias	VAV3; Guanine nucleotide exchange factor VAV3; VAV-3
Immunogen Species	Homo sapiens (Human)
Target Names	VAV3

Image



Western Blot analysis of 293T-UV cells using Vav3 Polyclonal Antibody

Usage

For Research Use Only. Not for use in diagnostic or therapeutic procedures.