




VTI1B Antibody

Product Code	CSB-PA006555
Storage	Upon receipt, store at -20°C or -80°C. Avoid repeated freeze.
Uniprot No.	Q9UEU0
Immunogen	Synthesized peptide derived from the Internal region of Human v-SNARE Vti1p.
Raised In	Rabbit
Species Reactivity	Human,Mouse,Rat
Tested Applications	WB, ELISA; WB:1:500-1:2000, ELISA:1:10000
Form	Liquid
Conjugate	Non-conjugated
Storage Buffer	Liquid in PBS containing 50% glycerol, 0.5% BSA and 0.02% sodium azide.
Purification Method	The antibody was affinity-purified from rabbit antiserum by affinity-chromatography using epitope-specific immunogen.
Isotype	IgG
Alias	VTI1B; VTI1; VTI1L; VTI1L1; VTI2; Vesicle transport through interaction with t-SNAREs homolog 1B; Vesicle transport v-SNARE protein Vti1-like 1; Vti1-rp1
Immunogen Species	Homo sapiens (Human)
Target Names	VTI1B

Image



Western Blot analysis of RAT-HEART cells using v-SNARE Vti1p Polyclonal Antibody

Usage

For Research Use Only. Not for use in diagnostic or therapeutic procedures.