




Phospho-CEBPA (T230) Antibody

Product Code	CSB-PA006574
Storage	Upon receipt, store at -20°C or -80°C. Avoid repeated freeze.
Uniprot No.	P49715
Immunogen	Synthesized peptide derived from Human C/EBP α around the phosphorylation site of T230.
Raised In	Rabbit
Species Reactivity	Human, Mouse, Rat
Tested Applications	WB, ELISA; WB:1:500-1:2000, ELISA:1:10000
Form	Liquid
Conjugate	Non-conjugated
Storage Buffer	Liquid in PBS containing 50% glycerol, 0.5% BSA and 0.02% sodium azide.
Purification Method	The antibody was affinity-purified from rabbit antiserum by affinity-chromatography using epitope-specific immunogen.
Isotype	IgG
Alias	CEBPA; CCAAT/enhancer-binding protein alpha; C/EBP alpha
Immunogen Species	Homo sapiens (Human)
Target Names	CEBPA

Image



Western Blot analysis of HepG2 cells using Phospho-C/EBP α (T230) Polyclonal Antibody