




Phospho-KAT5 (S86) Antibody

Product Code	CSB-PA006605
Storage	Upon receipt, store at -20°C or -80°C. Avoid repeated freeze.
Uniprot No.	Q92993
Immunogen	Synthesized peptide derived from Human TIP60 around the phosphorylation site of S86.
Raised In	Rabbit
Species Reactivity	Human, Mouse
Tested Applications	WB, ELISA; WB:1:500-1:2000, ELISA:1:40000
Form	Liquid
Conjugate	Non-conjugated
Storage Buffer	Liquid in PBS containing 50% glycerol, 0.5% BSA and 0.02% sodium azide.
Purification Method	The antibody was affinity-purified from rabbit antiserum by affinity-chromatography using epitope-specific immunogen.
Isotype	IgG
Alias	KAT5; HTATIP; TIP60; Histone acetyltransferase KAT5; 60 kDa Tat-interactive protein; Tip60; Histone acetyltransferase HTATIP; HIV-1 Tat interactive protein; Lysine acetyltransferase 5; cPLA(2)-interacting protein
Immunogen Species	Homo sapiens (Human)
Target Names	KAT5

Image



Western Blot analysis of Jurkat cells using Phospho-TIP60 (S86) Polyclonal Antibody

Usage For Research Use Only. Not for use in diagnostic or therapeutic procedures.