



Phospho-CENPA (S7) Antibody

Product Code	CSB-PA006641
Storage	Upon receipt, store at -20°C or -80°C. Avoid repeated freeze.
Uniprot No.	P49450
Immunogen	Synthesized peptide derived from Human CENP-A around the phosphorylation site of S7.
Raised In	Rabbit
Species Reactivity	Human
Tested Applications	IF, ELISA; IF:1:200-1:1000, ELISA:1:10000
Form	Liquid
Conjugate	Non-conjugated
Storage Buffer	Liquid in PBS containing 50% glycerol, 0.5% BSA and 0.02% sodium azide.
Purification Method	The antibody was affinity-purified from rabbit antiserum by affinity-chromatography using epitope-specific immunogen.
Isotype	IgG
Alias	CENPA; Histone H3-like centromeric protein A; Centromere autoantigen A; Centromere protein A; CENP-A
Immunogen Species	Homo sapiens (Human)
Target Names	CENPA
Usage	For Research Use Only. Not for use in diagnostic or therapeutic procedures.