



ERN2 Antibody

Product Code	CSB-PA006646
Storage	Upon receipt, store at -20°C or -80°C. Avoid repeated freeze.
Uniprot No.	Q76MJ5
Immunogen	Synthesized peptide derived from the C-terminal region of Human IRE1 α / β .
Raised In	Rabbit
Species Reactivity	Human, Mouse
Tested Applications	IHC, IF, ELISA; IHC:1:100-1:300, IF:1:200-1:1000, ELISA:1:10000
Form	Liquid
Conjugate	Non-conjugated
Storage Buffer	Liquid in PBS containing 50% glycerol, 0.5% BSA and 0.02% sodium azide.
Purification Method	The antibody was affinity-purified from rabbit antiserum by affinity-chromatography using epitope-specific immunogen.
Isotype	IgG
Alias	ERN2; IRE2; Serine/threonine-protein kinase/endoribonuclease IRE2; Endoplasmic reticulum-to-nucleus signaling 2; Inositol-requiring protein 2; hIRE2p; Ire1-beta; IRE1b
Immunogen Species	Homo sapiens (Human)
Target Names	ERN2
Usage	For Research Use Only. Not for use in diagnostic or therapeutic procedures.