




SH2B2 Antibody

Product Code	CSB-PA006660
Storage	Upon receipt, store at -20°C or -80°C. Avoid repeated freeze.
Uniprot No.	O14492
Immunogen	Synthesized peptide derived from the Internal region of Human APS.
Raised In	Rabbit
Species Reactivity	Human,Mouse,Rat
Tested Applications	WB, ELISA; WB:1:500-1:2000, ELISA:1:20000
Form	Liquid
Conjugate	Non-conjugated
Storage Buffer	Liquid in PBS containing 50% glycerol, 0.5% BSA and 0.02% sodium azide.
Purification Method	The antibody was affinity-purified from rabbit antiserum by affinity-chromatography using epitope-specific immunogen.
Isotype	IgG
Alias	SH2B2; APS; SH2B adapter protein 2; Adapter protein with pleckstrin homology and Src homology 2 domains; SH2 and PH domain-containing adapter protein APS
Immunogen Species	Homo sapiens (Human)
Target Names	SH2B2

Image



Western Blot analysis of HuvEc cells using APS Polyclonal Antibody

Usage For Research Use Only. Not for use in diagnostic or therapeutic procedures.