




Phospho-PDLIM5 (Y251) Antibody

Product Code	CSB-PA006664
Storage	Upon receipt, store at -20°C or -80°C. Avoid repeated freeze.
Uniprot No.	Q96HC4
Immunogen	Synthesized peptide derived from Human PDLIM5 around the phosphorylation site of Y251.
Raised In	Rabbit
Species Reactivity	Human
Tested Applications	WB, ELISA; WB:1:500-1:2000, ELISA:1:10000
Form	Liquid
Conjugate	Non-conjugated
Storage Buffer	Liquid in PBS containing 50% glycerol, 0.5% BSA and 0.02% sodium azide.
Purification Method	The antibody was affinity-purified from rabbit antiserum by affinity-chromatography using epitope-specific immunogen.
Isotype	IgG
Alias	PDLIM5; ENH; L9; PDZ and LIM domain protein 5; Enigma homolog; Enigma-like PDZ and LIM domains protein
Immunogen Species	Homo sapiens (Human)
Target Names	PDLIM5

Image



Western Blot analysis of K562 cells using Phospho-PDLIM5 (Y251) Polyclonal Antibody

Usage

For Research Use Only. Not for use in diagnostic or therapeutic procedures.