



Phospho-MYL12A/MYL12B/MYL9 (T17/S18) Antibody

Product Code	CSB-PA006671
Storage	Upon receipt, store at -20°C or -80°C. Avoid repeated freeze.
Uniprot No.	P19105/O14950/P24844
Immunogen	Synthesized peptide derived from Human MLC-2 around the phosphorylation site of T17/S18.
Raised In	Rabbit
Species Reactivity	Human, Mouse, Rat
Tested Applications	IHC, ELISA; IHC:1:100-1:300, ELISA:1:10000
Form	Liquid
Conjugate	Non-conjugated
Storage Buffer	Liquid in PBS containing 50% glycerol, 0.5% BSA and 0.02% sodium azide.
Purification Method	The antibody was affinity-purified from rabbit antiserum by affinity-chromatography using epitope-specific immunogen.
Isotype	IgG
Alias	MYL12A; MLCB; MRLC3; RLC; Myosin regulatory light chain 12A; MLC-2B; Myosin RLC; Myosin regulatory light chain 2; nonsarcomeric; Myosin regulatory light chain MRLC3; MYL12B; MRLC2; MYLC2B; Myosin regulatory light chain 12B; MLC-2A; MLC-2; M
Immunogen Species	Homo sapiens (Human)
Target Names	MYL12A/MYL12B/MYL9
Usage	For Research Use Only. Not for use in diagnostic or therapeutic procedures.