




IGF2BP2 Antibody

Product Code	CSB-PA006676
Storage	Upon receipt, store at -20°C or -80°C. Avoid repeated freeze.
Uniprot No.	Q9Y6M1
Immunogen	Synthesized peptide derived from the Internal region of Human IGF2-BP2.
Raised In	Rabbit
Species Reactivity	Human
Tested Applications	WB, ELISA; WB:1:500-1:2000, ELISA:1:20000
Form	Liquid
Conjugate	Non-conjugated
Storage Buffer	Liquid in PBS containing 50% glycerol, 0.5% BSA and 0.02% sodium azide.
Purification Method	The antibody was affinity-purified from rabbit antiserum by affinity-chromatography using epitope-specific immunogen.
Isotype	IgG
Alias	IGF2BP2; IMP2; VICKZ2; Insulin-like growth factor 2 mRNA-binding protein 2; IGF2 mRNA-binding protein 2; IMP-2; Hepatocellular carcinoma autoantigen p62; IGF-II mRNA-binding protein 2; VICKZ family member 2
Immunogen Species	Homo sapiens (Human)
Target Names	IGF2BP2

Image



Western Blot analysis of HuvEC cells using IGF2-BP2 Polyclonal Antibody

Usage For Research Use Only. Not for use in diagnostic or therapeutic procedures.