




CCT8 Antibody

Product Code	CSB-PA006683
Storage	Upon receipt, store at -20°C or -80°C. Avoid repeated freeze.
Uniprot No.	P50990
Immunogen	Synthesized peptide derived from the Internal region of Human TCP-1 θ .
Raised In	Rabbit
Species Reactivity	Human, Mouse
Tested Applications	WB, ELISA; WB:1:500-1:2000, ELISA:1:5000
Form	Liquid
Conjugate	Non-conjugated
Storage Buffer	Liquid in PBS containing 50% glycerol, 0.5% BSA and 0.02% sodium azide.
Purification Method	The antibody was affinity-purified from rabbit antiserum by affinity-chromatography using epitope-specific immunogen.
Isotype	IgG
Alias	CCT8; C21orf112; CCTQ; KIAA0002; T-complex protein 1 subunit theta; TCP-1-theta; CCT-theta; Renal carcinoma antigen NY-REN-15
Immunogen Species	Homo sapiens (Human)
Target Names	CCT8

Image



Western Blot analysis of HeLa cells using TCP-1 θ Polyclonal Antibody

Usage

For Research Use Only. Not for use in diagnostic or therapeutic procedures.