



Phospho-NFAT5 (S1197) Antibody

Product Code	CSB-PA006687
Storage	Upon receipt, store at -20°C or -80°C. Avoid repeated freeze.
Uniprot No.	O94916
Immunogen	Synthesized peptide derived from Human NFAT5 around the phosphorylation site of S1197.
Raised In	Rabbit
Species Reactivity	Human, Mouse
Tested Applications	IF, ELISA; IF:1:200-1:1000, ELISA:1:10000
Form	Liquid
Conjugate	Non-conjugated
Storage Buffer	Liquid in PBS containing 50% glycerol, 0.5% BSA and 0.02% sodium azide.
Purification Method	The antibody was affinity-purified from rabbit antiserum by affinity-chromatography using epitope-specific immunogen.
Isotype	IgG
Alias	NFAT5; KIAA0827; TONEBP; Nuclear factor of activated T-cells 5; NF-AT5; T-cell transcription factor NFAT5; Tonicity-responsive enhancer-binding protein; TonE-binding protein; TonEBP
Immunogen Species	Homo sapiens (Human)
Target Names	NFAT5
Usage	For Research Use Only. Not for use in diagnostic or therapeutic procedures.