




Phospho-PTGES3 (S113) Antibody

Product Code	CSB-PA006697
Storage	Upon receipt, store at -20°C or -80°C. Avoid repeated freeze.
Uniprot No.	Q15185
Immunogen	Synthesized peptide derived from Human p23 around the phosphorylation site of S113.
Raised In	Rabbit
Species Reactivity	Human,Mouse,Rat
Tested Applications	WB, IHC, ELISA; WB:1:500-1:2000, IHC:1:100-1:300, ELISA:1:10000
Form	Liquid
Conjugate	Non-conjugated
Storage Buffer	Liquid in PBS containing 50% glycerol, 0.5% BSA and 0.02% sodium azide.
Purification Method	The antibody was affinity-purified from rabbit antiserum by affinity-chromatography using epitope-specific immunogen.
Isotype	IgG
Alias	PTGES3; P23; TEBP; Prostaglandin E synthase 3; Cytosolic prostaglandin E2 synthase; cPGES; Hsp90 co-chaperone; Progesterone receptor complex p23; Telomerase-binding protein p23
Immunogen Species	Homo sapiens (Human)
Target Names	PTGES3

Image



Western Blot analysis of 3T3 cells using Phospho-p23 (S113) Polyclonal Antibody

Usage

For Research Use Only. Not for use in diagnostic or therapeutic procedures.