




CAPN9 Antibody

Product Code	CSB-PA006713
Storage	Upon receipt, store at -20°C or -80°C. Avoid repeated freeze.
Uniprot No.	O14815
Immunogen	Synthesized peptide derived from the Internal region of Human Calpain 9.
Raised In	Rabbit
Species Reactivity	Human
Tested Applications	WB, ELISA; WB:1:500-1:2000, ELISA:1:5000
Form	Liquid
Conjugate	Non-conjugated
Storage Buffer	Liquid in PBS containing 50% glycerol, 0.5% BSA and 0.02% sodium azide.
Purification Method	The antibody was affinity-purified from rabbit antiserum by affinity-chromatography using epitope-specific immunogen.
Isotype	IgG
Alias	CAPN9; NCL4; Calpain-9; Digestive tract-specific calpain; New calpain 4; nCL-4; Protein CG36
Immunogen Species	Homo sapiens (Human)
Target Names	CAPN9

Image



Western Blot analysis of HuvEc cells using Calpain 9 Polyclonal Antibody