




OR511 Antibody

Product Code	CSB-PA006717
Storage	Upon receipt, store at -20°C or -80°C. Avoid repeated freeze.
Uniprot No.	Q13606
Immunogen	Synthesized peptide derived from the Internal region of Human Olfactory receptor 511.
Raised In	Rabbit
Species Reactivity	Human
Tested Applications	WB, IF, ELISA; WB:1:500-1:2000, IF:1:200-1:1000, ELISA:1:10000
Form	Liquid
Conjugate	Non-conjugated
Storage Buffer	Liquid in PBS containing 50% glycerol, 0.5% BSA and 0.02% sodium azide.
Purification Method	The antibody was affinity-purified from rabbit antiserum by affinity-chromatography using epitope-specific immunogen.
Isotype	IgG
Alias	OR511; OLF1; Olfactory receptor 511; Olfactory receptor OR11-159; Olfactory receptor-like protein OLF1
Immunogen Species	Homo sapiens (Human)
Target Names	OR511

Image



Western Blot analysis of 293T cells using Olfactory receptor 511 Polyclonal Antibody

Usage

For Research Use Only. Not for use in diagnostic or therapeutic procedures.