




MRPS30 Antibody

Product Code	CSB-PA006740
Storage	Upon receipt, store at -20°C or -80°C. Avoid repeated freeze.
Uniprot No.	Q9NP92
Immunogen	Synthesized peptide derived from the C-terminal region of Human MRP-S30.
Raised In	Rabbit
Species Reactivity	Human
Tested Applications	WB, ELISA; WB:1:500-1:2000, ELISA:1:10000
Form	Liquid
Conjugate	Non-conjugated
Storage Buffer	Liquid in PBS containing 50% glycerol, 0.5% BSA and 0.02% sodium azide.
Purification Method	The antibody was affinity-purified from rabbit antiserum by affinity-chromatography using epitope-specific immunogen.
Isotype	IgG
Alias	MRPS30; PDCE9; BM-047; 28S ribosomal protein S30; mitochondrial; MRP-S30; S30mt; Programmed cell death protein 9
Immunogen Species	Homo sapiens (Human)
Target Names	MRPS30

Image



Western Blot analysis of LOVO cells using MRP-S30 Polyclonal Antibody