




MALT1 Antibody

Product Code	CSB-PA006741
Storage	Upon receipt, store at -20°C or -80°C. Avoid repeated freeze.
Uniprot No.	Q9UDY8
Immunogen	Synthesized peptide derived from the Internal region of Human MALT1.
Raised In	Rabbit
Species Reactivity	Human, Mouse
Tested Applications	WB, IHC, IF, ELISA; WB:1:500-1:2000, IHC:1:100-1:300, IF:1:200-1:1000, ELISA:1:10000
Form	Liquid
Conjugate	Non-conjugated
Storage Buffer	Liquid in PBS containing 50% glycerol, 0.5% BSA and 0.02% sodium azide.
Purification Method	The antibody was affinity-purified from rabbit antiserum by affinity-chromatography using epitope-specific immunogen.
Isotype	IgG
Alias	MALT1; MLT; Mucosa-associated lymphoid tissue lymphoma translocation protein 1; MALT lymphoma-associated translocation; Paracaspase
Immunogen Species	Homo sapiens (Human)
Target Names	MALT1

Image



Western Blot analysis of HUVEC cells using MALT1 Polyclonal Antibody

Usage For Research Use Only. Not for use in diagnostic or therapeutic procedures.