



KIF2C Antibody

Product Code	CSB-PA006755
Storage	Upon receipt, store at -20°C or -80°C. Avoid repeated freeze.
Uniprot No.	Q99661
Immunogen	Synthesized peptide derived from Human KIF2C around the non-phosphorylation site of S95.
Raised In	Rabbit
Species Reactivity	Human, Mouse
Tested Applications	IF, ELISA; IF:1:200-1:1000, ELISA:1:5000
Form	Liquid
Conjugate	Non-conjugated
Storage Buffer	Liquid in PBS containing 50% glycerol, 0.5% BSA and 0.02% sodium azide.
Purification Method	The antibody was affinity-purified from rabbit antiserum by affinity-chromatography using epitope-specific immunogen.
Isotype	IgG
Alias	KIF2C; KNSL6; Kinesin-like protein KIF2C; Kinesin-like protein 6; Mitotic centromere-associated kinesin; MCAK
Immunogen Species	Homo sapiens (Human)
Target Names	KIF2C
Usage	For Research Use Only. Not for use in diagnostic or therapeutic procedures.