




Phospho-CHEK1 (S317) Antibody

Product Code	CSB-PA006768
Storage	Upon receipt, store at -20°C or -80°C. Avoid repeated freeze.
Uniprot No.	O14757
Immunogen	Synthesized peptide derived from Human Chk1 around the phosphorylation site of S317.
Raised In	Rabbit
Species Reactivity	Human,Mouse,Rat
Tested Applications	WB, IHC, ELISA; WB:1:500-1:2000, IHC:1:100-1:300, ELISA:1:5000
Form	Liquid
Conjugate	Non-conjugated
Storage Buffer	Liquid in PBS containing 50% glycerol, 0.5% BSA and 0.02% sodium azide.
Purification Method	The antibody was affinity-purified from rabbit antiserum by affinity-chromatography using epitope-specific immunogen.
Isotype	IgG
Alias	CHEK1; CHK1; Serine/threonine-protein kinase Chk1; CHK1 checkpoint homolog; Cell cycle checkpoint kinase; Checkpoint kinase-1
Immunogen Species	Homo sapiens (Human)
Target Names	CHEK1

Image



Western Blot analysis of 293T cells using Phospho-Chk1 (S317) Polyclonal Antibody

Usage

For Research Use Only. Not for use in diagnostic or therapeutic procedures.