



CIT Antibody

Product Code	CSB-PA006777
Storage	Upon receipt, store at -20°C or -80°C. Avoid repeated freeze.
Uniprot No.	O14578
Immunogen	Synthesized peptide derived from the Internal region of Human Citron.
Raised In	Rabbit
Species Reactivity	Human, Mouse
Tested Applications	IHC, ELISA; IHC:1:100-1:300, ELISA:1:40000
Form	Liquid
Conjugate	Non-conjugated
Storage Buffer	Liquid in PBS containing 50% glycerol, 0.5% BSA and 0.02% sodium azide.
Purification Method	The antibody was affinity-purified from rabbit antiserum by affinity-chromatography using epitope-specific immunogen.
Isotype	IgG
Alias	CIT; CRK; KIAA0949; STK21; Citron Rho-interacting kinase; CRK; Serine/threonine-protein kinase 21
Immunogen Species	Homo sapiens (Human)
Target Names	CIT
Usage	For Research Use Only. Not for use in diagnostic or therapeutic procedures.