




CAPN10 Antibody

Product Code	CSB-PA006779
Storage	Upon receipt, store at -20°C or -80°C. Avoid repeated freeze.
Uniprot No.	Q9HC96
Immunogen	Synthesized peptide derived from the N-terminal region of Human Calpain 10.
Raised In	Rabbit
Species Reactivity	Human
Tested Applications	WB, ELISA; WB:1:500-1:2000, ELISA:1:40000
Form	Liquid
Conjugate	Non-conjugated
Storage Buffer	Liquid in PBS containing 50% glycerol, 0.5% BSA and 0.02% sodium azide.
Purification Method	The antibody was affinity-purified from rabbit antiserum by affinity-chromatography using epitope-specific immunogen.
Isotype	IgG
Alias	CAPN10; KIAA1845; Calpain-10; Calcium-activated neutral proteinase 10; CANP 10
Immunogen Species	Homo sapiens (Human)
Target Names	CAPN10

Image



Western Blot analysis of LOVO cells using Calpain 10 Polyclonal Antibody

Usage

For Research Use Only. Not for use in diagnostic or therapeutic procedures.