



ADCY6/ADCY5 Antibody

Product Code	CSB-PA006784
Storage	Upon receipt, store at -20°C or -80°C. Avoid repeated freeze.
Uniprot No.	O43306/O95622
Immunogen	Synthesized peptide derived from the C-terminal region of Human A Cyclase V/VI.
Raised In	Rabbit
Species Reactivity	Human, Mouse, Rat
Tested Applications	IHC, IF, ELISA; IHC:1:100-1:300, IF:1:200-1:1000, ELISA:1:20000
Form	Liquid
Conjugate	Non-conjugated
Storage Buffer	Liquid in PBS containing 50% glycerol, 0.5% BSA and 0.02% sodium azide.
Purification Method	The antibody was affinity-purified from rabbit antiserum by affinity-chromatography using epitope-specific immunogen.
Isotype	IgG
Alias	ADCY6; KIAA0422; Adenylate cyclase type 6; ATP pyrophosphate-lyase 6; Adenylate cyclase type VI; Adenylyl cyclase 6; Ca(2+)-inhibitable adenylyl cyclase; ADCY5; Adenylate cyclase type 5; ATP pyrophosphate-lyase 5; Adenylate cyclase type V;
Immunogen Species	Homo sapiens (Human)
Target Names	ADCY6/ADCY5
Usage	For Research Use Only. Not for use in diagnostic or therapeutic procedures.