



AKAP11 Antibody

Product Code	CSB-PA006794
Storage	Upon receipt, store at -20°C or -80°C. Avoid repeated freeze.
Uniprot No.	Q9UKA4
Immunogen	Synthesized peptide derived from the C-terminal region of Human AKAP 220.
Raised In	Rabbit
Species Reactivity	Human
Tested Applications	IHC, IF, ELISA; IHC:1:100-1:300, IF:1:200-1:1000, ELISA:1:10000
Form	Liquid
Conjugate	Non-conjugated
Storage Buffer	Liquid in PBS containing 50% glycerol, 0.5% BSA and 0.02% sodium azide.
Purification Method	The antibody was affinity-purified from rabbit antiserum by affinity-chromatography using epitope-specific immunogen.
Isotype	IgG
Alias	AKAP11; AKAP220; KIAA0629; A-kinase anchor protein 11; AKAP-11; A-kinase anchor protein 220 kDa; AKAP 220; hAKAP220; Protein kinase A-anchoring protein 11; PRKA11
Immunogen Species	Homo sapiens (Human)
Target Names	AKAP11
Usage	For Research Use Only. Not for use in diagnostic or therapeutic procedures.