




AKAP10 Antibody

Product Code	CSB-PA006795
Storage	Upon receipt, store at -20°C or -80°C. Avoid repeated freeze.
Uniprot No.	O43572
Immunogen	Synthesized peptide derived from the N-terminal region of Human AKAP 10.
Raised In	Rabbit
Species Reactivity	Human,Mouse,Rat
Tested Applications	WB, IHC, IF, ELISA; WB:1:500-1:2000, IHC:1:100-1:300, IF:1:200-1:1000, ELISA:1:10000
Form	Liquid
Conjugate	Non-conjugated
Storage Buffer	Liquid in PBS containing 50% glycerol, 0.5% BSA and 0.02% sodium azide.
Purification Method	The antibody was affinity-purified from rabbit antiserum by affinity-chromatography using epitope-specific immunogen.
Isotype	IgG
Alias	AKAP10; A-kinase anchor protein 10; mitochondrial; AKAP-10; Dual specificity A kinase-anchoring protein 2; D-AKAP-2; Protein kinase A-anchoring protein 10; PRKA10
Immunogen Species	Homo sapiens (Human)
Target Names	AKAP10

Image



Western Blot analysis of mouse cells using AKAP 10 Polyclonal Antibody

Usage For Research Use Only. Not for use in diagnostic or therapeutic procedures.