




DDX20 Antibody

Product Code	CSB-PA006796
Storage	Upon receipt, store at -20°C or -80°C. Avoid repeated freeze.
Uniprot No.	Q9UHI6
Immunogen	Synthesized peptide derived from the Internal region of Human Gemin3.
Raised In	Rabbit
Species Reactivity	Human, Mouse
Tested Applications	WB, ELISA; WB:1:500-1:2000, ELISA:1:40000
Form	Liquid
Conjugate	Non-conjugated
Storage Buffer	Liquid in PBS containing 50% glycerol, 0.5% BSA and 0.02% sodium azide.
Purification Method	The antibody was affinity-purified from rabbit antiserum by affinity-chromatography using epitope-specific immunogen.
Isotype	IgG
Alias	DDX20; DP103; GEMIN3; Probable ATP-dependent RNA helicase DDX20; Component of gems 3; DEAD box protein 20; DEAD box protein DP 103; Gemin-3
Immunogen Species	Homo sapiens (Human)
Target Names	DDX20

Image



Western Blot analysis of 293 cells using Gemin3 Polyclonal Antibody

Usage For Research Use Only. Not for use in diagnostic or therapeutic procedures.