




CA5B Antibody

Product Code	CSB-PA006799
Storage	Upon receipt, store at -20°C or -80°C. Avoid repeated freeze.
Uniprot No.	Q9Y2D0
Immunogen	Synthesized peptide derived from the C-terminal region of Human CA VB.
Raised In	Rabbit
Species Reactivity	Human,Mouse,Rat
Tested Applications	WB, IHC, IF, ELISA; WB:1:500-1:2000, IHC:1:100-1:300, IF:1:200-1:1000, ELISA:1:10000
Form	Liquid
Conjugate	Non-conjugated
Storage Buffer	Liquid in PBS containing 50% glycerol, 0.5% BSA and 0.02% sodium azide.
Purification Method	The antibody was affinity-purified from rabbit antiserum by affinity-chromatography using epitope-specific immunogen.
Isotype	IgG
Alias	CA5B; Carbonic anhydrase 5B; mitochondrial; Carbonate dehydratase VB; Carbonic anhydrase VB; CA-VB
Immunogen Species	Homo sapiens (Human)
Target Names	CA5B

Image



Western Blot analysis of 3T3 cells using CA VB Polyclonal Antibody

Usage

For Research Use Only. Not for use in diagnostic or therapeutic procedures.