




CHRM1 Antibody

Product Code	CSB-PA006806
Storage	Upon receipt, store at -20°C or -80°C. Avoid repeated freeze.
Uniprot No.	P11229
Immunogen	Synthesized peptide derived from the Internal region of Human mAChR M1.
Raised In	Rabbit
Species Reactivity	Human, Mouse, Rat
Tested Applications	WB, ELISA; WB:1:500-1:2000, ELISA:1:5000
Form	Liquid
Conjugate	Non-conjugated
Storage Buffer	Liquid in PBS containing 50% glycerol, 0.5% BSA and 0.02% sodium azide.
Purification Method	The antibody was affinity-purified from rabbit antiserum by affinity-chromatography using epitope-specific immunogen.
Isotype	IgG
Alias	CHRM1; Muscarinic acetylcholine receptor M1
Immunogen Species	Homo sapiens (Human)
Target Names	CHRM1

Image



Western Blot analysis of HepG2 cells using mAChR M1 Polyclonal Antibody

Usage

For Research Use Only. Not for use in diagnostic or therapeutic procedures.