




Phospho-TIRAP (Y86) Antibody

Product Code	CSB-PA006823
Storage	Upon receipt, store at -20°C or -80°C. Avoid repeated freeze.
Uniprot No.	P58753
Immunogen	Synthesized peptide derived from Human TIRAP around the phosphorylation site of Y86.
Raised In	Rabbit
Species Reactivity	Human, Mouse
Tested Applications	WB, IHC, ELISA; IHC:1:100-1:300, ELISA:1:5000
Form	Liquid
Conjugate	Non-conjugated
Storage Buffer	Liquid in PBS containing 50% glycerol, 0.5% BSA and 0.02% sodium azide.
Purification Method	The antibody was affinity-purified from rabbit antiserum by affinity-chromatography using epitope-specific immunogen.
Isotype	IgG
Alias	TIRAP; MAL; Toll/interleukin-1 receptor domain-containing adapter protein; TIR domain-containing adapter protein; Adaptor protein Wyatt; MyD88 adapter-like protein
Immunogen Species	Homo sapiens (Human)
Target Names	TIRAP

Image



Western blot analysis of K562 using p-TIRAP (Y86) antibody.

Usage For Research Use Only. Not for use in diagnostic or therapeutic procedures.