




ADCY9 Antibody

Product Code	CSB-PA006850
Storage	Upon receipt, store at -20°C or -80°C. Avoid repeated freeze.
Uniprot No.	O60503
Immunogen	Synthesized peptide derived from the N-terminal region of Human A Cyclase IX.
Raised In	Rabbit
Species Reactivity	Human
Tested Applications	WB, IHC, ELISA; WB:1:500-1:2000, IHC:1:100-1:300, ELISA:1:10000
Form	Liquid
Conjugate	Non-conjugated
Storage Buffer	Liquid in PBS containing 50% glycerol, 0.5% BSA and 0.02% sodium azide.
Purification Method	The antibody was affinity-purified from rabbit antiserum by affinity-chromatography using epitope-specific immunogen.
Isotype	IgG
Alias	ADCY9; KIAA0520; Adenylate cyclase type 9; ATP pyrophosphate-lyase 9; Adenylate cyclase type IX; Adenylyl cyclase 9
Immunogen Species	Homo sapiens (Human)
Target Names	ADCY9

Image



Western Blot analysis of HuvEc cells using A Cyclase IX Polyclonal Antibody

Usage

For Research Use Only. Not for use in diagnostic or therapeutic procedures.